

IASTM tool : Mustache

FASCIQ®

Convex edge rounded L

Technique

- Release layer/septa
- Deep tissue treatment
- Pin & stretch
- Effleurage

Region

- All body parts

Concave sharp edge L

Technique

- Scanning
- Superficial treatment

Region

- Arms
- Legs
- Shoulders
- Neck

Convex sharp edge XS

Technique

- Scanning
- Superficial

Region

- Foot
- Hand
- Achilles tendon
- Knee
- AC joint
- Elbow



Spot small

Technique

- Triggerpoint
- Release layers

Region

- All body parts
- Foot & hands

Concave round/sharp edge L

Technique

- Scanning
- Superficial
- Pin & stretch
- Effleurage

Region

- Lower leg
- Shoulder
- Arm

Convex sharp edge L

Technique

- Scanning
- Superficial treatment

Region

- All body parts



THYSOL